

# Label-free recognition of non-activated and activated human T cells by Quantative Phase Imaging

I.Andreyev<sup>1</sup>, J. Van Hauwermeiren<sup>2</sup>, J.Barbau<sup>2</sup>, L.Desmecht<sup>2</sup>, N.Szita<sup>1</sup>

1) Department of Biochemical Engineering, UCL, Gower Street, London, UK

2) Ovizio Cell Imaging, Rue du Bourbon 100 bus 2, 1180 Uccle, BELGIUM

E-mail: i.andreyev@ucl.ac.uk




## Abstract

There is a gap in current methods for evaluating quality and durability of cells in immunotherapeutic products like CAR-T. Ability to identify and characterize immune cells on single-cell level without labelling would be highly valuable for biomedical research and clinical trials where unmodified cells are required. We are developing analytical applications for label-free monitoring of human T cells activation and viral transduction by optical means using image based cell profiling methods. An application for detection of T cells activation by Quantative Phase Imaging (QPI) was developed. We showed that stimulated cells could be statistically distinguished from non-stimulated by optical morphology fingerprints obtained from their phase images. Miniaturization of cell bioprocessing platforms and performing cell Quality Control (QC) analysis in a fast and low-cost way would be a benefit. Microfluidic devices were tested and Quantative Phase Imaging was applied for detecting of non-activated and activated T cells in microfluidic chip. This presents a first step towards high throughput analytical platforms for CAR-T cells.


## Results

Label free Cell QC method was developed for monitoring of activation T cells by QPI




- We are developing analytical methods for label-free monitoring of activation and viral transduction of human T cells by optical means
- Non-activated and activated human T cells were imaged by qMod digital holographic camera (Ovizio Imaging Systems) on slides for 2D flow cytometry (Chemometrics) and shown on Fig. A
- Optical morphology fingerprints were obtained from cell holographic images for each cell being imaged. Then by Principal Component Analysis a few of these fingerprints, which have the most difference for non-activated and activated cells, were identified for further use
- A plot of Phase Optical Volume (POV), which is directly proportional to cell dry mass, versus cell perimeter is shown on a Fig. B above
- It was shown that stimulated T cells could be statistically distinguished from non-stimulated cells by optical morphology fingerprints obtained from their phase images. Standard statistical approach for evaluation of any new method for detection of a disease by biomarkers was used
- Evaluation of this method using Cell Phase Optical Volume as a parameter is shown on the insert to the Fig. C. Evaluation of sensitivity/specificity of detection of activated vs non-activated T cells taken from a second donor showed similar statistical output

Label free recognition of non-activated and activated T cells on a chip by QPI




- The diffusion microfluidic cell culture chip was designed and manufactured
- It is optically transparent and built from Glass and Si
- Melted agar gel is loaded into the chip and guided by hydrophobic Au strips forming vertical gel membranes
- Gel membranes diffuse gases and nutrients in cell chamber
- It was shown that activated and non-activated T cells could be distinguished by QPI and the monitoring of activation of immunotherapeutic T cells *on a chip* is feasible




Figures up – schematic and prototype of the cell culture chip used for imaging; Figures down/left – typical phase images of non-activated and activated T cells; Figures down/right – plot of cell morphology fingerprints (cell Phase OV vs Perimeter) obtained from phase images, on the insert – statistics calculated from Phase Optical Volume distributions for non-activated and activated T cells

Can we in the future monitor viral loading in human T cells?

- Infected cells show changes on cell phase images (shown for infected Insected cells)



- Viral loading monitoring of infected Insected cells



- Allowing optimization of viral procedures
- Can be done for any virus produced in cells

- Changes obtained on phase images of infected Insected cells are believed might be obtained for a viral transduced CAR-T cells allowing monitoring of viral loading in immunotherapeutic T cells

## Future work

- Developing Critical Analytics applications for in-line QC of therapeutic cells using Quantative PI
- Developing relevant methods and techniques of high throughput image/spectra based cell profiling with emphasis on Deep Learning and Big Data
- Apply Raman spectroscopy for phenotyping of primary T/CAR-T cells
- Developing microfluidic cell culture for therapeutic T cells processing and analysis

## Acknowledgements

Funding from the UK Engineering & Physical Sciences Research Council (EPSRC) for the Future Targeted Healthcare Manufacturing Hub hosted at University College London with UK university partners is gratefully acknowledged (Grant Reference: EP/P006485/1) Financial and in-kind support from the consortium of industrial users and sector organisations is also acknowledged

In collaboration with

