



iLine F PRO

LIVE CELL ANALYZER

Optimize your bioprocess efficiency
through real-time cell monitoring



VISIT [OVIZIO.COM](https://www.ovizio.com)



Shed light on your cells

The Ovizio iLine F PRO provides you with improved insight into your cell-based processes by bringing **real-time monitoring of Cell Quality Attributes (CQA's)**.

Based on **label-free analysis**, our analyzer ensures continuous monitoring of your cells status' without the manual handling usually required.

Secured connection to your bioreactor is ensured by our BioConnect, a **closed-loop cartridge** system that reduces the risk of contamination and further reduces manual handling.


Cell detection is made possible through **machine-learning algorithms**, enabling continuous tracking of multiple parameters such as viability, bead count, cell activation profile, cell infection status, etc.

Our **user-friendly software** integrates remote-monitoring and automation features for improved cell-monitoring.

Our solution has followed numerous regulatory controls for **cGMP compliance readiness**, making it the ideal candidate for bioprocessing applications in cell & gene therapy applications.


The result:

- Real-time information on cell culture status
- Drastically reduced hands-on time
- Maximized availability of biological sample
- Accelerated development of cell-based process



Features

Benefits



Automated analyzer
Closed-loop system

Achieve continuous cell monitoring
Drastically reduce handling operations
Drastically reduce risk of external contamination
Eliminate wasted biological material



Machine Learning
algorithms

Achieve reproducible measurements
Identify multiple Cell Quality Attributes (CQAs)
Customize your cell classification



Remote monitoring

Continuous cell monitoring
Early detection of cell-growth issues
Alert triggering options



Open platform
communication

Optimize your bioprocesses & feeding strategy
Reduce time-to-harvest


iLine F PRO

The iLine F PRO enables **label-free cell analysis** through holographic imaging. The analyzer is easy to use and only features one button for on/off operations. It is equipped with a built-in computer and a 7 inch touchscreen display to enable continuous display of the ongoing cell monitoring.

Data can easily be backed up through a built-in USB port or sent to a dedicated server through the integrated ethernet port. An integrated QR code reader enables BioConnect identification for improved traceability.




BioConnect



Our BioConnect enables **closed-loop connection** between the iLine F PRO analyzer and the bioreactor. Through Luer lock or welding, it can connect with all common types of bioreactors, such as rocking-motion bags, re-usable and disposable benchtop bioreactors or stainless-steel bioreactors. Its assembly meets many ISO and FDA norms such as ISO 10 993, ISO 14 644, CFR 21 210 and CFR 21 211. It has been designed to be cGMP compliant ready for manufacturing vaccines or cell- & gene-therapies.

Machine Learning Algorithm

Our machine learning algorithm will scan each image and identify all visible objects, such as cells and beads. For each of these objects, it will measure multiple features and combine them into a biological signature.





This signature will be utilized for cell classification and enables identification of different cell culture parameters such as:



- Cell viability
- Cell count
- Magnetic bead count
- Morphological changes



Our software uses machine learning to identify object features and classify cells accordingly

	Situation A	Situation B
		
Object Features		
Diameter	10	7
Circularity	10	8
Membrane Regularity	10	1
...
Biological Classifier :	(10,10,0,...)	(7,8,1,...)

OsOne Software



Example of CHO Cells cultivated in Dasgib bioreactor for 7 days.

The OsOne software is at the center of Ovizio's analyzers. Designed to deliver an attractive user experience, OsOne is built for easy data acquisition and thorough quantitative data analysis of suspension cell cultures. The intuitive interface and wizard enable a rapid understanding and navigation through the app. Improved user traceability has been integrated to ensure compliance with FDA norm CFR21/11.

Specifications

Ovizio reader

Type	iLine F holographic transmission microscope
Light source	Partially coherent monochromatic LED 630 nm (red)
Sensor	C-MOS camera 2456 (H) x 2058 (V) – 8 bits
Analyzer objective	Olympus LUCPLFLN20X analyzer objective 20x – 0.45 NA
Total magnification	22.2X
Horizontal resolution	1.5 μm
Field of view	318 μm x 318 μm
Acquisition rate	15 fps max
Data volume	128 Mb per measurement point. Note: volume may vary depending on selected cell line
Input Power	100–240 V AC ~ 47–63 Hz – 45 W Typical, 90 W Max

Integration capabilities

Device	400 (15.7" – L) x 223 (9.2" – W) x 423 (16.7" – H) mm – 21 kg
Shipment box	515 x 320 x 550 mm

OsOne specifications

OPC ready, automation via RESTful web services
Compliant with CFR21/11 regulations

Data output: object count, viability, viable cell density, total cell density...

BioConnect

Cell concentration range	Cell counting range of 0.5 to 10 Mcells/mL The device provides a \pm 10% cell density variation
Typical flow rate	From 4 to 12 mL/min
Cell travel time	5 to 25 min. Varies according to specific setup (pump profile, total ...)
Tubing length	180 cm of silicon 1 x 3 mm + 2 * 40 cm of PVC 3.1 x 4.2 mm
Welding interface	PVC tubing (3.1 x 4.2 mm ID x OD)
Addressable bioreactor volume ¹	From 300 mL to 10.0 L total process
Retention volume	4.8 mL to 8.3 mL Note: Retention volume of BioConnect only, other tubing added between vessel and BioConnect will increase these values
Sterilization	Gamma irradiation (25 kGy) and individual packaging
Flow cell channel height	Between 150 and 200 μm
Hygrometry	20–80% – non-condensing
Temperature	15–25 °C
Shelf life	6 months at storage temperature from 15 to 30 °C
Bioreactor inner pressure limitation	Max 7.25 [500] psig [mbarg]
Regulatory compliance	All fluid-contact materials have passed current USP Biological Test or ISO 10993 certification

Note: these specifications may change without notice.

¹Tested range for efficient pumping and monitoring

